



35. PTEV SERVER ROOM

A room shall be reserved for the computer servers of the PTEV Therapy Control System (TCS) equipment. All cables from the PTEV data network shall emerge from under the raised floor in this room for connection to the PTEV server through adequate patch panel (refer to Fig.67.01). The following is required:

1. Located at the treatment level.
2. As close as possible to the MCR.
3. Minimum of 12 m² (130 ft²).
4. Raised floor (antistatic computer floor), free height 30 cm (1'), flush with treatment level floor.
5. Suspended dismountable ceiling (height 2.5 m (8'-3") measured from finished floor).
6. Independent and stand alone air conditioning unit. This computer room shall be equipped with an air-cleaning system adapted for a proper server room. See figure 61.01 for heat dissipation.
7. Dust shall not penetrate in a quantity to interfere with satisfactory operation of the apparatus or to impair safety.
8. The D/BT shall provide:
 - 8.1. 1 desk (150 cmx80 cm, 5'x2'-7") and 2 chairs
 - 8.2. 1 safe box and cupboard for backup tapes (±40 cm x 80 cm , 1'-4"x2'-8")
 - 8.3. 1 switch cabinet with min 12 U free for PTEV. A full height cabinet 80 cm (2'-8") depth is recommended.
9. The D/BT shall provide a high grade epoxy paint finish. The exposed concrete floor slab below the raised computer flooring shall receive an industrial type epoxy coating and the walls below the raised computer flooring shall be epoxy painted to control dust. The epoxy paint type shall be submitted to PETV for approval prior purchasing.
10. Lighting shall be on a dimming system to prevent computer screen glare.
11. The door will be lockable with auto-closing and a possibility to keep it open. The access keys (main and spare) will be handed over to PTEV only.
12. No sprinkler or wet pipes, drains in the ceiling or raised floor.
13. For network details refer to Chapter VI : Other technics Set : Network and for UPS details Figure 52.32



REV:	DATE:	MODIFICATION:	DRAFTSMAN:	CHECKED BY:	VALIDATED BY:
A	30/04/15	Original issue	LCHEN	QBA	PV

MATERIAL: --
 SCALE: 1/25 (A3)
 DIMENSIONS: mm
 TOLERANCES: --

PROJECT: PROTON THERAPY
 SPROJECT: TATA HBTF MUMBA

II.
 ROOMS

Other rooms

TITLE:
 PTEV Server Room

07.42.33.

35.01 A

THIS DOCUMENT IS THE PROPERTY OF IBA AND SHALL NOT BE ALTERED, COPIED OR COMMUNICATED WITHOUT THE WRITTEN CONSENT OF IBA.
 IBA
 S.A. ION BEAM APPLICATIONS

36. MACHINE SHOP

If the scope of the project includes a Machine Shop, it shall be equipped with the necessary equipment to produce treatment patient specific devices such as compensators. In addition to this, the customer may want to add equipment in this room for maintenance purposes. Size will be determined in accordance with equipment to be located in this room. The following is required:

1. If required a hatch in the ceiling of the machine shop, refer to *Figure 50.01-1, Table 2: Definition of Hatches*, will be used as access for repairing or installing beam line equipment. Or an equivalent access through the building.
 2. Double doors, giving access to the treatment corridor (minimum doors aperture: 185 cm (6') width x 240 cm (8') height. If the milling machine must be transported through these doors, refer to the dimensions given hereunder.
 3. Compressed air outlets at each milling machine + 1 additional.
 4. Industrial size double sink with cold and warm water.
 5. Cold water outlet at each milling machine.
 6. An electrical connection box and the connection for (depending contract):
 - 6.1 Each Milling Machine (400V, +/- 30KVA for each milling machine).
 - 6.2 Lathe
 - 6.3 Band Saw
 - 6.4 Additional Customer Equipment
 7. Floor drain at each Milling machine.
 8. Noise insulation/sound proofing in walls and ceiling.
 9. Industrial heavy duty shelving will be provided by the D/BT. A detailed drawing of the shelves will be made by PTEV after the dimensions of the room are fixed.
 10. Fire proof cabinet for oil and flammable liquids will be provided by the D/BT.
 11. Special care need to be taken to ensure access into the building for the equipment.

Example milling machine (Current size may be subject to small change):

Weight: 4800 Kg (10500 lbs)
 Width: 270 cm (8'-10")
 Depth: 320 cm (10'-6")
 Height: 270 cm (8'-10")
- Allow a minimum of 60 cm (2') clearances all around each milling machine for maintenance purpose and above for lighting and cabling

REV:	DATE:	MODIFICATION:	DRAFTSMAN:	CHECKED BY:	VALIDATED BY:
1	30/04/15	Original issue	LCHEN	QBA	PV

MATERIAL:	--
SCALE:	(A3)
DIMENSIONS:	mm
TOLERANCES:	--

PROJECT: PROTON THERAPY
 SPROJECT: TATA HBTF MUMBA

II. ROOMS

Other rooms


TITLE:
 Machine Shop

07.42.33.

36.01 A

THIS DOCUMENT IS THE PROPERTY OF IBA AND SHALL NOT BE ALTERED WITHOUT THE WRITTEN CONSENT OF IBA.

S.A. ION BEAM APPLICATIONS



37. OTHER PTEV ROOMS

1. MAINTENANCE ROOM

A room shall be reserved to make some maintenance on PTE small electrical and mechanical parts. The following is required:

- 1.1. Minimum 20 m² (215 ft²).
- 1.2. Connected or close to the Machine Shop or the PTEV office.
- 1.3. A sink with cold and hot water.
- 1.4. A desk and a chair.
- 1.5. Compressed air outlets.
- 1.6. Outlets 230 V (110V US).
- 1.7. Built-in caseworks and cabinets
- 1.8. Ductwork and piping in the ceiling areas should be kept to a minimum in order to maximize the space.
- 1.9. Fire proof cabinet for oil and flammable liquids will be provided by the D/BT.
- 1.10. The doors will be lockable with auto-closing and a possibility to keep them open.

2. PTEV SPARE PARTS STORAGE ROOM

Two rooms shall be reserved for the storage of PTE spare parts, accessible and equipped from the BOD.

One room at the treatment floor level

- Minimum of 50 m² (560 ft²)

One room at the upper level (PSR level)

- Minimum of 50 m² (560 ft²)
- As close as possible to the PSR

The following is required:

- 2.1. Dust shall not penetrate in a quantity to interfere with satisfactory operation of the apparatus or to impair safety. This must be operational since BOD.
- 2.2. Industrial heavy duty shelving and drawers (Vidmar type).
- 2.3. The final drawing of the shelves and drawers will be made by D/BT after the dimensions of the room are fixed and provided to PTEV at the latest 3 months prior to BOD.
- 2.4. Pallet jack capability (+/- 300 kg (600 lbs))
- 2.5. Standard double lockable door access with auto-closing and a possibility to keep them open.
- 2.6. 1 desk and chair in each room.
- 2.7. Lighting will be positioned according to shelves position.
- 2.8. Easy access to the service elevator with double doors all the way of routing.
- 2.9. Ceiling and walls will have a light dustproof painting and floors will have an industrial epoxy finish.

3. PTEV OFFICES AND TEMPORARY SPACES

Offices shall be reserved for the PTEV engineers and operators.

During construction:


- 3 months before BOD: Office with furniture for 3 people.
- 2 months before BOD: Storage for survey tools (see also ML-3)
- 2 weeks before BOD:
- Office with furniture for 8 to 10 people until move to permanent offices.
- Office with furniture for 3 people.

From BOD and after acceptance, permanent offices for maintenance (operation phases):

- One office space for 10 - 12 people with furniture (lockable and ergonomic desks, chairs and cabinets for 14 people) - Minimum of 85 m² (920 ft²).
- Two offices for 3 people (min 11m² (120 ft²) each) equipped with standard furniture - Located at the treatment level, as close as possible to the MCR.
- A storage and copy room (5 m² (55 ft²)) with network plugs.
- A temporary conference and meeting room from the beginning of the installation period up until the end of the acceptance phase. Equipped with desks and chairs for 12 people

4. ELEVATOR

At least one elevator must be available to transport a load of minimum 2500 kg (5500 Lbs) between different levels of the building. This shall be used occasionally for the transportation of spare parts. The elevator entry door shall be a minimum clear width of 140cm(4'-6") x245cm (8') high. The cab size would be at minimum 215cm (7') deep by 185cm (6') wide.



REV:	A	DATE:	30/04/15	MODIFICATION:	Original Issue
CHECKED BY:	PV	VALIDATED BY:			
CHECKED BY:	OBA				
DRAFTSMAN:	LCHEN				

MATERIAL:	-
SCALE:	(A3)
DIMENSIONS:	mm
TOLERANCES:	-

PROJECT: PROTON THERAPY
 SPROJECT: TATA HBTF MUMBAI

II. ROOMS

Other rooms

TITLE:
 Other PTEV rooms

07.42.33.

37.01 A

III. INSTALLATION

50. INSTALLATION

ACRONYMS

- BOD: Building occupancy date
- ESS: Energy selection system
- D/BT: Design/building team
- PTEV: Proton therapy equipment vendor
- BTS: Beam transport system
- PTE: Proton Therapy Equipment
- Ri: Route i
- PPS: Patient Positioning System


TABLE OF CONTENTS:

- Routes definition 50.01-1
- Routes definition 50.01-2
- Alignment sighting 50.02

REFER ALSO TO:

- Fig. 01.11 (Beam Transport Components & Isocenters)
- Fig. 01.11-2, Table 2: Cyclotron Main Components
- Fig. 01.11-2, Table 3: ESS/BTS Main Components
- Fig. 02.11, Table 1: Gantry Structure Main Components & Nozzle by room
- Fig. 02.13 (Gantry Loads)
- Fig. 02.14 (PPS Loads)
- Fig. 02.51 (Gantry Treatment Room Dimensions-Installation Beam & Lifting Points)
- Fig. 02.52 (Gantry Nozzle Anchor for Lifting Device Installation)
- Fig. 02.53 (PPS Pit Embeds)
- Fig. 32.01, Table 2: Cooling Room Main Components
- Fig. 33.01, Table 1: MCR Main Components
- Fig. 34.01, Table 1: Treatment Control Room Equipment
- Fig. 52.11-2 (Cabinets Numbering)
- Fig. 90.01 (Building Completion Requirements)
- Fig. 90.02 (Milestones)

THIS DOCUMENT IS THE PROPERTY OF IBA AND SHALL NOT BE ALTERED COPIED OR COMMUNICATED WITHOUT THE WRITTEN CONSENT OF IBA.



REV: A	DATE: 30/04/15	MODIFICATION: Original issue	DRAFTSMAN: LCHEN	CHECKED BY: QBA	VALIDATED BY: PV

MATERIAL: -
 SCALE: (A3)
 DIMENSIONS: mm
 TOLERANCES: -

PROJECT: PROTON THERAPY
 PROJECT: TATA HBTF MUMBAI

III. INSTALLATION

TITLE: Presentation

07.42.33.

50.00 A

1. BENCHMARKS TO BE PROVIDED BY THE D/BT

The PTEV will layout the equipment alignment grid and delineate it with permanent monument markers to be installed after the BOD. The D/BT shall provide a primary benchmark (Foundation Target Insert (FTI) as supplied by Hubbs Machine & Manufacturing, Inc) at a location surveyed with respect to the building grid system. The PTEV will define the equipment alignment grid including elevation with respect to this primary benchmark. The primary benchmark will be located in the floor near the corner between the North and West walls of the BTS.

This primary benchmark shall be documented by the D/BT. The PTEV shall also install secondary benchmarks that are surveyed with respect to the above primary benchmark as indicated in *Fig.01.11-1 - Beam Transport Components & Isocenters*

The exact coordinates of the secondary benchmarks will be determined during the equipment installation.

All benchmarks need to remain accessible all the time without putting any equipment in front of them. In case of a finish wall need to be erected, access panels will be foreseen to keep access to the benchmarks.

2. ALIGNMENT SIGHTING PROVISIONS

The PTEV will require the D/BT to make provision for alignment sightings through the shield walls.

The quantities, size, and location of those sighting holes are specified in *Fig.50.02*. When several alignment view axes intercept, it must be assured that a common volume is visible along all these view axes. All equipment or building components that may block the view must be removable such that a clear line of sight can be achieved during equipment realignment.

The D/BT must provide manually removable shielding plugs for the sighting holes (this can be accomplished using interlocking removable concrete shielding blocks). The design will be validated by PTEV and D/BT shielding consultant.

Closing of sighting holes by D/BT after PTEV agreement.

3. ALIGNMENT OF EQUIPMENT WITH RESPECT TO THE BUILDING

3.1 THERMAL STABILIZATION REQUIREMENTS

The Facility HVAC systems in all PTE equipment areas shall be operational prior to Building Occupation Date (BOD) by PTEV, per phased BOD dates. Please also refer to *Fig 90.02*.



Before starting alignment of the equipment the Facility must be stabilized to room temperature and all structural walls and slabs must have reached thermal equilibrium. Therefore, the start date for the alignment process is linked to closing the building temporary hatches and completing a thermal soak period with the HVAC system in normal operation.

3.2 FLOOR REQUIREMENTS

The beam line equipment extending from the cyclotron through to the Gantries must be aligned to a level horizontal plane with high precision. The first step in accomplishing this is using shims or equivalent under the beam line stand pads to level the stands and bring them to the correct height.

The second step is to fill the gap between the pads and floor with high strength, non-shrink, non-metallic grout (5000 PSI (35 MPa) in 3 days). See also the building completion requirements in *Fig.90.01*.

The grouting shall be done by D/BT on PTEV request, and sealed and epoxy painted.



Filling with concrete cylinders by D/BT. Proposal shall be submitted to PTEV for approval

REV:	A	DATE:	30/04/15
MODIFICATION:	Original Issue	DRAFTSMAN:	LCHEN
CHECKED BY:	QBA	VALIDATED BY:	PV

MATERIAL: —
 SCALE: 1/250 (A3)
 DIMENSIONS: mm
 TOLERANCES: —

PROJECT: PROTON THERAPY
 SPROJECT: TATA HBTf MUMBAI

III. INSTALLATION

TITLE:
 Alignment Sighting

07.42.33.

50.02 A

IV. COOLING / GASES

51. COOLING & GASES

ACRONYMS

BOD:	Building Occupancy Date
D/BT:	Design Building team
PTEV:	Proton Therapy Equipment Vendor
BTS:	Beam Transport System
PSR:	Power Supply Room
WCR:	Water Cooling Room
PTE:	Proton Therapy Equipment
CC:	Cooling Circuit
GTR:	Gantry Treatment Room

TABLE OF CONTENTS:

- Introduction	51.01
- Cooling Schematic Distribution in the Building	51.02
- Cooling Schematic: Tables	51.03
- General Cooling Requirements	51.11
- Cooling Ducts in the Building	51.12-1/2
- Water Cooling Room & PSR Water Distribution	51.13
- Main Coil water conditioner #1	51.21
- Cyclotron vacuum pump #2	51.22
- Cyclo, ESS, BTS & Gantry water conditioner #3	51.23
- Compressed air, dry nitrogen & gases requirements	51.31

REFER ALSO TO:

Figure 00.31-2 (Facility General Layout: Basement & Upper Level)

Figure 32.01 (Water Cooling Room)

Refer also to **Chapter II (ROOMS)** for architectural specifications of the rooms

THIS DOCUMENT IS THE PROPERTY OF IBA AND SHALL NOT BE ALTERED COPIED OR COMMUNICATED WITHOUT THE WRITTEN CONSENT OF IBA.



S.A. ION BEAM APPLICATIONS

REV:	A	DATE:	30/04/15	MODIFICATION:	Original Issue
DRAFTSMAN:	LCHEN	CHECKED BY:	OBA	VALIDATED BY:	PV

MATERIAL: -
 SCALE: (A3)
 DIMENSIONS: mm
 TOLERANCES: -

PROJECT: PROTON THERAPY
 SPROJECT: TATA HBTF MUMBAI

IV. COOLING & GASES

TITLE:
Presentation

07.42.33.

51.00 A

INTRODUCTION

This chapter is related to water cooling system and gases necessary for PTE. First the cooling water distribution is detailed in schematic of the cooling circuits and in the tables. Then integration in the building and water conditioners details. Finally, we'll see the requirements related to the gases at the end of the chapter.

The water distribution is based on 3 independent circuits:

- Cyclotron Main Coils cooling circuit (#1)
- Cyclotron Vacuum Pumps cooling circuit (#2)
- Cyclotron (remaining)/ESS/BTS/GTR/PSR cooling circuit (#3)

Figure 51.12-2 (Cooling Ducts in the building) gives the implantation of the cooling circuits in the building. An overview of the flow rates necessary for cooling the different parts of the equipment is given in Table 1: Hydraulic Characteristics figure 51.03.

As a convention, the Primary circuit is the PTEV side and the Secondary circuit is the D/BT side.

The 3 related water conditioners are located in the Water Cooling room. They will be provided and installed by the PTEV. See also Figure 00.31-2 (Facility General Layout: Basement & Upper Level) and Figure 32.01 (Water Cooling Room) for implantation in the building.

CIRCUIT DESIGN


FOR EACH COOLING CIRCUIT, THE FOLLOWING IS REQUIRED:

1. Refer to Figure 51.02 and 51.03 for pressure, temperature, flow and other calculation inputs.
2. The secondary circuit of the water-conditioner is a closed loop that feeds from the building chilled water supply. The secondary loop shall satisfy the parameters in Figure 51.03.
3. The D/BT shall provide and locate the primary and secondary water hookups in order to allow an easy connection between these hookups and the interface flanges of the water conditioner. The location and types of the flanges are shown in Figure 51.21, Figure 51.22 and Figure 52.23 and in embedded Table 2: Flanges Specification for the Water Conditioners. These flanges will be supplied to the D/BT by the PTEV.
4. The final piping sections going from the hookups to the PTEV water conditioners must be completed by the D/BT, after the equipment is installed in the water cooling room.
5. The secondary cooling circuit shall be run with a bypass before connecting it to the PTEV water conditioner in order to remove any construction waste in the pipes that could damage the heat exchanger. A pipeline strainer is needed with valves on both sides.
6. See embedded Table 1 for water quality.
7. The piping and any part in contact with water must be compatible with high purity de-ionized water. Materials compatibility will be checked by D/BT.

TABLE1: Water Quality			
S.No.	Parameter	Standards	Rationale/Remarks
1.	pH range	6.5-8.5	The range is conducive for propagation of aquatic species and restoring natural system
2.	Turbidity	30 NTU	Reasonably clear water for Recreation, Aesthetic appreciation and Industrial cooling purposes.
3.	Dissolved Iron	0.5 mg/l or less	It is desirable to have the (as Fe) collective concentration of dissolved Fe and Mn less or equal to 3 x 10-5 pound/foot3 to avoid scaling effect.
4.	Dissolved Manganese(as Mn)	3 x 10-5 pound/foot3	
5.	water conductivity	< 1 µS/cm (< 2.54 µS/inch)	

Table 2: Flanges Specification for the Water Conditioners

Item	Value for DN 50 (mm/inch)	Value for DN 100 (mm/inch)	Value for DN 150 (mm/inch)
Diameter A	165 / 6.50	220 / 8.66	285 / 11.2
Diameter B	125 / 4.92	180 / 7.09	240 / 9.4
Diameter C	~55 / 2.17	~110 / 4.33	~160 / 6.3
Diameter D	4 x 18 / 0.71	8 x 18 / 0.71	8 x 22 / 0.87



REV: A	DATE: 30/04/15	MODIFICATION: Original Issue	DRAFTSMAN: LCHEN	CHECKED BY: QBA	VALIDATED BY: PV
--------	----------------	------------------------------	------------------	-----------------	------------------

MATERIAL: -
 SCALE: (A3)
 DIMENSIONS: mm
 TOLERANCES: -

PROJECT: PROTON THERAPY
 SPROJECT: TATA HBTF MUMBAI

IV. COOLING & GASES

TITLE: Introduction

07.42.33.

51.01 A




REV:	A	DATE:	30/04/15	MODIFICATION:	Original Issue
DRAFTSMAN:	LCHEN	CHECKED BY:	OBA	VALIDATED BY:	PV

MATERIAL: (A3)
SCALE: mm
DIMENSIONS: mm
TOLERANCES: -

Legend

- Compressed air
 - Supply Water
 - Return Water
 - Reduction Valve
 - Release valve, 8 bars
 - Compressed air connection
- Extra tube in the concrete installed as back-up if initial pipes are clogged or leaking
- Sediment Filter, 5 µm
 - 1" Y Strainer, 250 µm, stainless steel mesh, PN10 bars
 - 1.5" Y Strainer, 250 µm, stainless steel mesh, PN10 bars
 - 2" Y Strainer, 250 µm, stainless steel mesh, PN10 bars
 - 2" Y Strainer, 250 µm, stainless steel mesh, PN20 bars
 - 1" Y Strainer, 500 µm, stainless steel mesh, PN10 bars
 - 1" Y Strainer, 500 µm, stainless steel mesh, PN20 bars
 - 2" Y Strainer, 500 µm, stainless steel mesh, PN10 bars
 - 2.5" Y Strainer, 250 µm, stainless steel mesh, PN10 bars
 - Y Strainer, 250 µm, stainless steel mesh, PN10 bars, Size to be determined by D/BT.
 - 1" Y Strainer, 1 mm, stainless steel mesh, PN10 bars
 - 1" Y Strainer, 1 mm, stainless steel mesh, PN20 bars
 - 2" Y Strainer, 1 mm, stainless steel mesh, PN20 bars
 - 4" Y Strainer, 1 mm, stainless steel mesh, PN10 bars
- (1) Filter and its isolating valve are temporary equipments that can be removed when the return piping is very cleaned of any impurities.
 - (2) Any water in excess must be evacuated out in conformity with the regional law.
 - (3) Some filter type may not be applicable for this configuration.



Exact position of the connections see Fig.51.12-2



Table 1: Hydraulic characteristics

	Cydo / ESS / BTS / Gantry / PSR	Main coil	Cydo vacuum system
Input differential pressure	Bar	4 to 6	
Supply Temperature	°C	5° - 35°	
Flow	m³/h	> 2 m³/h	
	Secondary circuit (Chiller side)		
	Interface Circuit 3	Interface Circuit 1	Interface Circuit 2
Water conditioner max static pressure	4	4	4
Input differential pressure	1	1	0.6
Supply Temperature	5 < t° < 15	5 < t° < 15	5 < t° < 15
Flow	110	20	ODP: 0.6
Max. Power to be removed	750	210	ODP: 5
	Primary circuit (ba Equipment side)		
	Circuit 3.1	CC1CYC, Circuit 1.1: Main Coil,	CC2CYC, Circuit 2
Piping test static pressure	12	25	6
Water conditioner max static pressure	8	18	4
Max flow	133	8	ODP: 0.6
Max water return temperature	40	50	35
Building Pipe diameter	8	2	1
	CC1PSR, Circuit 3.2, Power supply Room	CC2PSR, Circuit 1.2: Main Coil Power Supply	
Piping test static pressure	12	35(*)	
Water conditioner max static pressure	8	6	
Max flow	13	1.5	
Max water return temperature	40	35	
Building Pipe diameter	3(**)	1	
	(*) This value is related to case of dysfunction of the pressure regulator.		
	(**) This value is related to case of dysfunction of the pressure regulator.		
	CC3CYC, Circuit 3.3, Cycle Miscellaneous		
Piping test static pressure	12		
Water conditioner max static pressure	8		
Max flow	20		
Max water return temperature	40		
Building Pipe diameter	2		
	Circuit 3.4, ESS/BTS/Gantry		
Piping test static pressure	12		
Water conditioner max static pressure	8		
Max flow	600		
Max water return temperature	40		
Building Pipe diameter	8		
	Circuit 3.4, Sub-circuits		
	See figure 51.02		

Table 2: Water Distribution Legend

Connection ID	Equipment	Fitting size	PN	Fitting type	Fitting material	Pipes material
CCMKUP		1/2"	10	BSP Fem + ECM-RK	Brass	PVC-C (chlorinated polyvinyl chloride, high pressure) pipes is acceptable, nevertheless stainless steel is recommended
CC1SEC	Main Coil Water Conditioner	2"	10	Flange DN50	Brass	
CC2SEC	Cydo vacuum conditioner	1"	10	BSP Fem + ECM-RK	Brass	
CC3SEC	Cydo/Ess/Gantry / PSR	4"	10	Flange DN100	Brass	
CC11PRIM	Main Coil	2"	20	Flange DN50	Stainless steel	Stainless steel
CC12PRIM	Main Coil PS	1"	20	BSP Fem + ECM-RK	Stainless steel	Stainless steel
CC2PRIM	Cydo Vacuum pump	1"	10	BSP Fem + ECM-RK	Brass	PVC-C (chlorinated polyvinyl chloride, high pressure) pipes is acceptable, nevertheless stainless steel is recommended
CC31PRIM	Cyclotron misc. (FPA + manifolds)	Supply: 4" Return: 6"	10	Supply: Flange DN100 Return: Flange DN150	Brass	
CC1CYC	Cyclotron main coils	2"	20	Flange DN50	Stainless steel	Stainless steel
CC2CYC	Cyclotron vacuum pumps	1"	10	BSP thread, female	Brass	PVC-C (chlorinated polyvinyl chloride, high pressure) pipes is acceptable, nevertheless stainless steel is recommended
CC3CYC	Cyclotron misc. (FPA + manifolds)	2"	10	BSP thread, female	Brass	
CC1PSR	Power supplies in PSR	3"	10	BSP thread, female	Brass	
CC2PSR	Main Coil Power supply in PSR	1"	20	BSP thread, female	Stainless steel	Stainless steel
CC3ESS	ESS	1.5"	10	BSP thread, female	Brass	
CC1SLI	Static line TR1 (B1...1) + 2 quads of static line for next TR	1.5"	10	BSP thread, female	Brass	PVC-C (chlorinated polyvinyl chloride, high pressure) pipes is acceptable, nevertheless stainless steel is recommended
CC1TR1	GTR1	2.5"	10	BSP thread, female	Brass	
CC1TR2	GTR2	2.5"	10	BSP thread, female	Brass	
CC1TR3	GTR3	2.5"	10	BSP thread, female	Brass	

MATERIAL: -
SCALE: (A3)
DIMENSIONS: mm
TOLERANCES: -

PROJECT: PROTON THERAPY
SPROJECT: TATA HBTF MUMBAI

IV.
COOLING &
GASES

TITLE:
Cooling
Schematic: Tables

07.42.33.

51.03 A

REV: A	DATE: 30/04/15	MODIFICATION: Original Issue	DRAFTSMAN: LCHEN	CHECKED BY: CBA	VALIDATED BY: PV
--------	----------------	------------------------------	------------------	-----------------	------------------

Integration in the building

GENERAL REQUIREMENTS

- The D/BT is responsible for providing and installing all the pipes described in this section.
- The pipe sizes are estimated by PTEV according to the values defined and schematically represented in Figure 51.02 (Cooling Schematic Distribution in the Building). Nevertheless, the pipe sizing and the final configuration remain the responsibility of the D/BT. Any change shall be submitted to the PTEV for approval.
- The pipes shall be designed for maximum pressures drop of 0.5 bar (7.25 PSI) in each line at nominal flow rate.
- All the pipes shall be leak tight; the leak test will be performed by D/BT at the pressure referenced as piping test static pressure in Figure 51.02 (Cooling Schematic Distribution in the Building).
- The pipes shall be closed at their ends to avoid dust, concrete, obstructing the ducts until final connection.
- The water trenches will be for the PTEV use only. No other equipment can be installed in the trenches without the PTEV agreement. In that case, trenches should be made larger and the design adapted and validated by PTEV
- D/BT shall provide removable walk-on aluminum or galvanized steel metal gratings on all water trenches, similar to those requested for the wiring trenches.
- The trenches and lowest points in all rooms dedicated to PTEV shall be provided with drains connected to the main floor drain system.
- Metallic water conduits shall be grounded according to the applicable regulations.
- Pipe routing must be designed to minimize the radiation leakage.
- A purge outlet shall be provided at the lowest point of each circuit.
- When water pipes are embedded in the concrete, 2 spare pipes (supply/return) shall be added for back-up connections.

FLOOR DRAIN REQUIREMENTS

- Floor drains and/or sump capacity shall be provided in any area of the building entered by water pipes.
- The drains shall be sufficient to capture, without flooding, the total water volume estimated in Table 1: Estimated Total Water Capacity of Water Circuits in Cubic Meters (Gallons) for any cooling water circuit passing through that area.
- A sump filter shall also be installed on all the drains. In both trenches (along the beam line), drains shall be foreseen every 5 meters (±17').
- Floor drains shall receive same finish as the trenches and slabs.

BSP THREAD DEFINITION

Male thread fittings shall be tapered and female thread fittings shall be parallel. Threads according to the British Standard Institute - BS 21 "Specification for pipe threads for tubes and fittings where pressure-tight joints are made on the threads".

PARALLEL (STRAIGHT) THREADS TYPE:

BS parallel threads are also referred to as British Gas, British Pipe Parallel or Parallel Fastening Thread. Common symbols used for the thread: BSP, BSPP, BSSPI, BSPF, BSPG, PS, Rp, G. The BSPP (parallel) male will mate with a BSPP (parallel) female or a female port. Usually marked: G size ISO

TAPER THREADS TYPE:

BSP taper threads are also referred to as British Standard Taper Pipe, Pipe Taper or Conical Thread. Common symbols used for the thread: BSPT, BSPT_r, PT, KR, R, Rc Usually marked: Rc size ISO

Characteristics of BSPT taper threads:

- tapered thread 1°47'
- truncation of roots and crests are rounded
- 55° thread angle
- pitch is measured in threads per inch
- taper is 3/4" per foot or 1 in 16 on the diameter

Example of supplier for NPT to BSP adapter fittings: Mac Master-Carr. <http://www.mcmaster.com> (example: 1-1/2" NPT male x 1-1/2" BSPT male brass fitting ref. 5832T118)


THREAD DIMENSIONS: see Table 2

Table 1: Estimated Total Water Capacity of Water Circuits in Cubic Meters (Gallons) –to be updated

	Circuit #1 Conditioner 0.5 (132)	Circuit #2 Conditioner 0.03 (8)	Circuit #3 Conditioner 1.7 (450)
Conditioner	0.5 (132)	0.03 (8)	1.7 (450)
Piping	0.2 (53)	0.1 (26)	0.4 (106)
Equipment	0.45 (116)	0.23 (61)	0.4 (106)
Contingency	0.3 (80)	0.2 (53)	1 (254)
TOTAL	1.45 (383)	0.73 (193)	3.5 (925)

Table 2: THREAD DIMENSIONS

Pipe Size (inch)	Pitch (threads per inch)	Major Diameter OD (inch)
1/8"	28	0.383
1/4"	19	0.518
3/8"	19	0.656
1/2"	14	0.825
3/4"	14	1.041
7/8"	14	1.189
1"	11	1.309
1-1/4"	11	1.650
1-1/2"	11	1.882
2"	11	2.347
2-1/2"	11	2.900
3"	11	3.460



REV.	DATE:	MODIFICATION:	DRAFTSMAN:	CHECKED BY:	VALIDATED BY:
A	30/04/15	Original Issue	LCHEN	OBA	PV

MATERIAL: -
 SCALE: (A3)
 DIMENSIONS: mm
 TOLERANCES: -

PROJECT: PROTON THERAPY
 SPROJECT: TATA HBTF MUMBAI


IV. COOLING & GASES

TITLE:
 General Cooling Requirements

07.42.33.

Refer to figure 51.12-2

1. The main supply and return pipes of cooling circuit # 3 (as well as the compressed air service) shall enter in the BTS area at the point W located in this figure.
2. From this point, all the pipes will run in a trench called "Water Service Trench", located on the right hand side of the beam line when looking in the beam direction. It means that several pipes will run in this trench.
3. The supply and return pipes of cooling circuit #1, #2 and CC3CYC of circuit #3 will enter the cyclotron vault area through the ceiling shielding.




- Compressed air
- Nitrogen
- Oxygen
- Hydrogen

O2 and H2 outlets shall be installed as indicated on this drawing whether the bottles are installed in the pit or located elsewhere in the building


Caution: water from cooling system might be considered activated depending on local regulations. In this case, it shall be stored for decay during 1 day or according to local regulation.

Pipes going under the BTS shall be installed at minimum 300 mm (1') under treatment level, to prevent from interactions with table anchoring.



THIS DOCUMENT IS THE PROPERTY OF IBA AND SHALL NOT BE ALTERED OR COMMUNICATED WITHOUT THE WRITTEN CONSENT OF IBA.

SA. ION BEAM APPLICATIONS



REV:	A	DATE:	30/04/15	MODIFICATION:	Original Issue
DRAFTSMAN:	LCHEN	CHECKED BY:	QBA	VALIDATED BY:	PV

MATERIAL: —
 SCALE: 1/200 (A3)
 DIMENSIONS: mm
 TOLERANCES: —

PROJECT: PROTON THERAPY
 SPROJECT: TATA HBTF MUMBAI


IV. COOLING & GASES

TITLE:
 Cooling Ducts in the Building

07.42.33.

POWER SUPPLY ROOM WATER DISTRIBUTION


- 4 Pipes (Supply & Return) shall pass through the walls separating the water-cooling room and the power supply room. If the water pipes are embedded in the concrete, then 4 additional spare pipes shall be foreseen. They will emerge under the access floor and stop at the wall.
- The water characteristics for Cyclotron/ESS/BTS/PSR/Gantry distribution are also valid for the power supply room water distribution. The diameter of the pipes shall be according to Figure 51.02.



○ compressed air

THIS DOCUMENT IS THE PROPERTY OF IBA AND SHALL NOT BE ALTERED OR COMMUNICATED WITHOUT THE WRITTEN CONSENT OF IBA.

S.A. ION BEAM APPLICATIONS



REV:	A	DATE:	30/04/15	MODIFICATION:	Original Issue
DRAFTSMAN:	LCHEN	CHECKED BY:	QBA	VALIDATED BY:	PV

MATERIAL: —
 SCALE: 1/100 (A3)
 DIMENSIONS: mm
 TOLERANCES: —

PROJECT: PROTON THERAPY
 SPROJECT: TATA HBTF MUMBAI

IV. COOLING & GASES

TITLE:
 Water Cooling Room & Power Supply Room Water Distribution

07.42.33.

51.13 A

THIS DOCUMENT IS THE PROPERTY OF IBA AND SHALL NOT BE ALTERED, COPIED OR COMMUNICATED WITHOUT THE WRITTEN CONSENT OF IBA.



S.A. ION BEAM APPLICATIONS

REV:	A	DATE:	30/04/15	MODIFICATION:	Original issue
DRAFTSMAN:	LCHEN	CHECKED BY:	QBA	VALIDATED BY:	PV

MATERIAL: -
SCALE: 1/25 (A3)
DIMENSIONS: mm
TOLERANCES: -


PROJECT: PROTON THERAPY
SPROJECT: TATA HBTF MUMBAI

IV.
COOLING &
GASES


TITLE:
Main coil water
conditioner #1

07.42.33.


51.21 A



FRONT VIEW



SIDE VIEW



Top View

PORT	PURPOSE
S2	OUTLET SECONDARY, connected to Chiller
R2	INLET SECONDARY, connected to Chiller
R1	INLET PRIMARY, connected to Cyclotron
S1	OUTLET PRIMARY, connected to Cyclotron
S0	Fill in (MakeUp)
R0	Drain
S11	Outlet to IBA equipment in PSR
R11	Inlet from IBA equipment in PSR

Water conditioner weight:
• empty : 640 kg / 1410 lb
• loaded (water): 1050 kg / 2315 lb

INDICATIVE VALUE ONLY
SUBJECT TO SMALL MODIFICATIONS

REV:	A	DATE:	30/04/15	MODIFICATION:	Original Issue
DRAFTSMAN:	LCHEN	CHECKED BY:	OBA	VALIDATED BY:	PV

MATERIAL: (A3)
 SCALE: mm
 DIMENSIONS:
 TOLERANCES:






PROJECT: PROTON THERAPY
 SPROJECT: TATA HBTF MUMBAI

IV.
 COOLING &
 GASES

TITLE:
 Cyclotron vacuum
 pump #2

07.42.33.

51.22 A



PORT	FUNCTION	WATER FLOWS DIRECTION
S0	Fill in (MakeUp)	
R0	Drain	
S1	Outlet to IBA equipment	Water conditioner to IBA equipment
R1	Inlet from IBA equipment	IBA equipment to water conditioner
S2	Inlet from customer	Chiller to water conditioner
R2	Outlet to customer	Water conditioner to chiller

Water conditioner weight:
 • empty : 70 kg / 154 lb
 • loaded (water): 100 kg / 220 lb