

| | | | | | |
|------------|-------|-------------|----------|---------------|----------------|
| REV: | A | DATE: | 30/04/15 | MODIFICATION: | Original issue |
| DRAFTSMAN: | LCHEN | CHECKED BY: | QBA | VALIDATED BY: | PV |

MATERIAL: —
 SCALE: (A3)
 DIMENSIONS: mm
 TOLERANCES: —

PROJECT: PROTON THERAPY
 SPROJECT: TATA HBTF MUMBAI

**IV.
 COOLING &
 GASES**



TITLE:
 Cyclo, ESS, BTS
 & Gantry water
 conditioner #3

07.42.33.

51.23 A

| PORT | PURPOSE |
|------|--|
| R2 | OUTLET SECONDARY, connected to Chiller |
| S2 | INLET SECONDARY, connected to Chiller |
| R1 | INLET PRIMARY, connected to Cyclotron & BTS |
| S1 | OUTLET PRIMARY, connected to Cyclotron & BTS |
| R0 | Drain |
| S0 | Fill in (MakeUp) |

Water conditioner weight:
 • empty : 2600 kg / 5732 lb
 • loaded (water): 4900 kg / 10803 lb



COMPRESSED AIR

Compressed air is required to operate pneumatic actuators (for valves, beam stops), to vent parts of the PTE which are under vacuum and, finally, for the Gantry pneumatic brakes. The required compressed air outlets are indicated in Figure 51.12-1/2 (cooling ducts in the building).

General requirements:

- Compressed air is required in all the rooms where PTEV equipment is installed as indicated, except in the main control room and treatment control rooms. For the Power Supply Room and the Water Cooling Room the compressed air feed point shall be in a central location of the room. The D/BT shall provide a valve and threaded connection as described below in paragraph No. 9 at this location.
- Two independent and redundant compressors provided by the D/BT are required to ensure the availability of the Facility.
- At least one reservoir must be provided by the D/BT in order to have compressed air in case of a power failure. Air accumulator capacity in case of power supply failure to be globally 200 liters (52 gallons).
- Flow is discontinuous with peak flow 4 l/second (one gallon/second), but not more than once a minute, and normally less than once an hour.
- Required pressure is between 6 & 10 bars (87 psi & 145 psi).
- The D/BT will provide pressure test reports of the pressure vessel(s).
- Compressed air must be dust-free, oil free and industrial dry (per standard commercial practice for shop air) with a dew point of -20°C (-4°F).
- D/BT will provide all the outlets fittings in brass, 1" BSP female thread.
- D/BT will provide one independent circuit for the PTE side.

ESS/IBTS COMPRESSED AIR SPECIAL REQUIREMENTS

- A compressed air pipe provided by the D/BT (one inch diameter) will run in the cooling water service trench all along the beam line as shown in Figure 51.12-2 (Cooling Ducts & in the building) & Figure 51.02 (Cooling Schematic Distribution in the Building). This pipe will be used for purging the equipment of the beam line.

GANTRY COMPRESSED AIR SPECIAL REQUIREMENTS

- Clean and industrial dry air is required in the PPS pit and at the rear support roller assembly (as shown in Figure 51.12-1 (Cooling Ducts & in the building)).
- Dual airlines are required in order to enhance Facility availability.
- The average airflow rate per Gantry based on braking once every 1.5 minutes is 1 SCFM (1.7Nm³/h).
- Air to the junction box is to allow for pneumatic actuators on the Gantry. Valve and threaded connection shall be provided as described above.

MAINTENANCE ROOM

- If this room is part of the project scope, clean and industrial dry air outlets are required in the Maintenance Room.

DRY NITROGEN

Dry nitrogen outlet provided by the D/BT is required in order to vent the cyclotron PTE: see Figure 51.12-1/2 (Cooling Ducts & in the building). This outlet shall be connected to interchangeable nitrogen bottles (equipped with pressure regulator) to be placed by D/BT in a location dedicated for gas bottle storage, as required by the local regulation. PTEV recommends locating the nitrogen bottles in a rated closet.

General requirements:

- Industrial nitrogen 99.9 %
- b50 x 4 (4 bottles, a manifold and retention rack following local regulation)
- Pressure: between 1 & 2 bars (14.5 & 29 psi)
- Estimated consumption: 4 m³ (140 ft³) (STP) maximum per venting.
- Flow : > 2m³ in 5 minutes (70 ft³ in 5 min)
- Piping : Brass, copper, or stainless steel with label
- D/BT will provide the outlets fitting, 3/8" BSP female thread. Final connection to the cyclotron will be made by PTEV.

OXYGEN

Oxygen is required for the deflector inside the cyclotron. See Figure 51.12-1/2 (Cooling Ducts in the building) for outlet position.

If allowed by local regulation, the bottles will be located close to the cyclotron, at the degrader height, on the wall holding the manifolds. If not allowed by local regulation, the bottles will be placed as close as possible to the cyclotron vault (piping < 30m if possible).

General requirements:

- Purity: Oxygen N45 (i.e. 99.995 %) provided and installed by D/BT
- Bottle type : b20 x 2 (1 bottle connected and 1 spare, a manifold and retention rack following local regulation)
- Pressure regulation: between 1 & 3 bars (15 & 45 psi)
- Manometer on high and low pressure levels
- Flow : between 1 and 2 cm³/min (STP) (between 0.06 and 0.12 in³/min)
- Piping : Stainless Steel 1/16" with label
- Located in the pit with 300 cm ($\pm 10'$) free length to connect directly on the cyclotron.
- All equipment on the line will be in Stainless steel

HYDROGEN

Hydrogen is required inside of the Cyclotron. See Figure 51.12-1/2 (Cooling Ducts in the building) for outlet position.

If allowed by local regulation, the bottles will be located inside the cyclotron pit and there will be no interference with the building. If not allowed by local regulation, the bottles will be placed as close as possible to the cyclotron vault (piping < 30 m ($\pm 100'$) if possible).

General requirements:

- Purity: Hydrogen N60 (i.e. 99.9999 %)
- Bottle type : b10 x 2 (1 bottle connected and 1 spare, a manifold and retention rack following local regulation)
- Pressure regulation: between 1 & 3 bars (15 & 45 psi)
- Manometer on high and low pressure levels
- Flow : between 1 and 2 cm³ /min (STP) (between 0.06 and 0.12 in³/min)
- Piping : chromatography grade, steam cleaned stainless steel 1/16" with label on it
- Located in the pit with 5 m ($\pm 17'$) free length to connect directly on the cyclotron.
- All equipment on the line will be in Stainless steel
- D/BT shall ensure that pipes can be replaced easily in case of leakage or breakage

HELIUM

Helium is not used by the cyclotron. See Figure 51.12-1/2 (Cooling Ducts in the building) for outlet position.

General requirements:

- Purity: Helium N45 (i.e. 99.995 %)
- Bottle type : b5 or b10 x 2 (1 bottle on a carrier to roll along the whole beam line and 1 spare, a manifold and retention rack following local regulation)
- Provided with a flexible pipe of min. 4 m ($\pm 12'$) and a blow gun.
- Pressure regulation: between 1 & 1.5 bars (15 & 22 psi)
- A retention rack for 1 spare bottle, following local regulation. Provided and installed by D/BT in a place following local regulation.

| | | | | | |
|------------|-------|-------------|----------|---------------|----------------|
| REV.: | A | DATE: | 30/04/15 | MODIFICATION: | Original issue |
| DRAFTSMAN: | LOHEN | CHECKED BY: | OBA | VALIDATED BY: | PV |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |

| | |
|-------------|------|
| MATERIAL: | --- |
| SCALE: | (A3) |
| DIMENSIONS: | mm |
| TOLERANCES: | --- |

PROJECT: PROTON THERAPY
 SPROJCT:TATA HBTF MUMBAI

IV. COOLING & GASES

TITLE:

Compressed air,
dry nitrogen &
Gases
requirements

07.42.33.

51.31 A

V. ELECTRICAL SET

52. ELECTRICAL SET

ACRONYMS

| | |
|-------|---------------------------------|
| UPS: | Uninterruptible Power Supply |
| BTS: | Beam Transport System |
| PTE: | Proton Therapy Equipment |
| D/BT: | Design & Building Team |
| PSR: | Power Supply Room |
| PTEV: | Proton Therapy Equipment Vendor |
| WCR: | Water Cooling Room |
| CBCT: | Cone Beam Computed Tomography |
| PT: | Proton Therapy |
| ESS: | Energy Selection System |
| TSS: | Therapy Safety System |
| TCR: | Treatment Control Room |
| MCR: | Main Control Room |
| HK: | Hook Up |
| CT: | Cable Tray |
| PPS: | Patient Positioning System |
| BTB: | Building Terminal Box |
| SC: | Switching Cabinet |

TABLE OF CONTENTS:

| | |
|---|-----------|
| - Cabinets numbering | 52.11-1/2 |
| - Main power distribution schematic | 52.21 |
| - Interfaces requirements | 52.22 |
| - Interfaces examples | 52.23 |
| - Required electrical Hookups table | 52.24 |
| - Electrical Hookups locations | 52.31 |
| - UPS connections location | 52.32 |
| - Grounding connections location | 52.33 |
| - Cable trays & ducts requirements | 52.41 |
| - General comments and legend | 52.42 |
| - Schematic : cable distribution in the building | 52.43 |
| - Location of cable trays & duct ends: treatment level | 52.44 |
| - Location of cable trays & duct ends: Upper level & basement | 52.45 |
| - Cable trays & ducts views and details | 52.46 |
| - MCR cable trays | 52.47 |
| - Cable trays in cyclo vault, ESS/BTS | 52.48 |
| - PSR cabinets stands | 52.51 |

For reference documents, refer to figure 00.00

Refer also to **Chapter II (ROOMS)** for architectural specifications of the rooms



| | | | | | |
|------------|-------|-------------|----------|---------------|----------------|
| REV: | A | DATE: | 30/04/15 | MODIFICATION: | Original issue |
| DRAFTSMAN: | LCHEH | CHECKED BY: | QBA | VALIDATED BY: | PV |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |

| | |
|-------------|------|
| MATERIAL: | - |
| SCALE: | (A3) |
| DIMENSIONS: | mm |
| TOLERANCES: | - |

PROJECT: PROTON THERAPY
 PROJECT: TATA HBTF MUMBAI

V. ELECTRICAL SET

TITLE:
Presentation

07.42.33.

52.00 A

THIS DOCUMENT IS THE PROPERTY OF IBA AND SHALL NOT BE ALTERED, COPIED OR COMMUNICATED WITHOUT THE WRITTEN CONSENT OF IBA.



S.A. ION BEAM APPLICATIONS

| | | | | | | | | | | | |
|------|---|-------|----------|---------------|----------------|------------|-------|-------------|-----|---------------|----|
| REV: | A | DATE: | 30/04/15 | MODIFICATION: | Original Issue | DRAFTSMAN: | LCHEN | CHECKED BY: | QBA | VALIDATED BY: | PV |
|------|---|-------|----------|---------------|----------------|------------|-------|-------------|-----|---------------|----|


MATERIAL: —
 SCALE: 1/200 (A3)
 DIMENSIONS: mm
 TOLERANCES: —

PROJECT: PROTON THERAPY
 SPROJECT: TATA HBTF MUMBAI


V.
 ELECTRICAL SET

TITLE:
 Cabinets numbering

07.42.33.



UPPER LEVEL



TREATMENT LEVEL

THIS DOCUMENT IS THE PROPERTY OF IBA AND SHALL NOT BE ALTERED, COPIED OR COMMUNICATED WITHOUT THE WRITTEN CONSENT OF IBA.



| REV. | DATE: | MODIFICATION: | DRAFTSMAN: | CHECKED BY: | VALIDATED BY: |
|------|----------|----------------|------------|-------------|---------------|
| A | 30/04/15 | Original Issue | LCHEH | QBA | PV |

MATERIAL: ---
SCALE: 1/200 (A3)
DIMENSIONS: mm
TOLERANCES: ---

PROJECT: PROTON THERAPY
SPROJECT: TATA HBTF MUMBAI


V.
ELECTRICAL SET

TITLE: Cabinets numbering


07.42.33.

52.11-2 A 2/2

| Area Code | Cabinet Name | Location | W mm | D mm | H mm | Quantity |
|-----------|---|--------------------|------|------|------|----------|
| B001 | Cyclotron Main Cool PS | PSR | 1818 | 2050 | 1800 | 2000 |
| B017 | RFA Airside PS | PSR | 3150 | 2250 | 2100 | 2520 |
| B050 | ISide Heater | PSR | 600 | 230 | 600 | TRD |
| B118 | Power Distribution Cycle | PSR | 1800 | 420 | F | Stair |
| B091 | Power Distribution PT | PSR | 2100 | 420 | F | Stair |
| B220 | Quadrupoles General Transformer PS | PSR | 2100 | 800 | 2000 | 2700 |
| B320 | Dipoles General Transformer PS | PSR | 2100 | 800 | 2000 | 2700 |
| B430 | SAPS (Scanning Magnetic Power Supply) | PSR | 1800 | 500 | 1000 | Stair |
| B440 | SAPS EU Cab | PSR | 600 | 600 | 1100 | Stair |
| B455 | SAPS Switching | PSR | 1300 | 600 | 1000 | Stair |
| B470 | PSR Water Supply | CYC ROOM | 2100 | 800 | 1800 | 300 |
| A012 | Final Power Amplifier (100 KW) | PSR | 4600 | 600 | 1000 | 420 |
| B555 | IS5 IS5 1424 Switching System PS 1 | ESS | 600 | 800 | 1600 | 450 |
| B341 | Water/Vacuum/TE Electronic Cabinet | PSR | 600 | 800 | 1600 | 450 |
| D002 | BLU | PSR | 600 | 800 | 1600 | 450 |
| B500 | Tri-Power Supply Cabinet | PSR | 600 | 800 | 1600 | 450 |
| E001 | ES5 Electronic Cabinet (CPU) | ESS | 600 | 800 | 1600 | 450 |
| E002 | ES5 Electronic Cabinet (PT) | ESS | 600 | 800 | 1600 | 450 |
| E301 | BTS Electronic Cabinet 1 | BTS Vault | 425 | 220 | 600 | 200 |
| E302 | BTS Electronic Cabinet 2 | BTS Vault | 425 | 220 | 600 | 200 |
| E303 | BTS Electronic Cabinet 3 | BTS Vault | 425 | 220 | 600 | 200 |
| A036 | AP Distribution Cabinet | CYC ROOM | 600 | 210 | 800 | 100 |
| C022 | Powering Control Unit N1 | Pulling Room 1 | 600 | 800 | 1800 | 450 |
| B323 | X-Ray Generator Cabinet TR1 | Pulling Room 1 | 500 | 468 | 1101 | 348 |
| G308 | Gantry Driver Cabinet TR1 | Gantry Ph 1 | 800 | 420 | 1000 | F |
| G313 | Robot Controller Unit TR1 (PPS) | Pulling Room 1 | 810 | 600 | 1250 | F |
| G117 | Robot Controller Unit TR1 (PPS) HDC Box | Pulling Room 1 | 600 | 800 | 1800 | 450 |
| C021 | Treatment Control Unit N1 | Gantry 1 | 400 | 420 | 200 | 15 |
| G303 | Distribution Box TR1 | Gantry 1 | 600 | 800 | 1800 | 450 |
| G307 | Nozzle 1 Electronic Units Cabinet N1 | Gantry 1 | 553 | 780 | 767 | 200 |
| G308 | Nozzle 2 Electronic Units Cabinet N1 (OCE/UA) | Gantry 1 | 553 | 780 | 767 | 200 |
| G309 | Gantry 1 Electronic Cabinet N1 (BPS/BEU) | Gantry 1 | 553 | 780 | 767 | 200 |
| C011 | SICU1 | Pulling Room 1 | 1000 | 320 | 1000 | 70 |
| G306 | Gantry Derivatives Box TR1 - Param Plane1 | Gantry Ph 1 | 1200 | 420 | 1200 | 300 |
| G308 | Gantry Derivatives Box TR1 - Control Plane1 | Gantry Ph 1 | 1200 | 420 | 1200 | 300 |
| D023 | Powering Control Unit N2 | Pulling Room 2 | 600 | 800 | 1800 | 450 |
| R321 | X-Ray Generator Cabinet TR2 | Gantry Ph 2 | 500 | 468 | 1101 | 348 |
| G306 | Gantry Driver Cabinet TR2 | Gantry Ph 2 | 800 | 420 | 1000 | F |
| G313 | Robot Controller Unit TR2 (PPS) | Pulling Room 2 | 810 | 600 | 1250 | F |
| S017 | Robot Controller Unit TR2 (PPS) HDC Box | Pulling Room 2 | 600 | 800 | 1800 | 450 |
| D021 | Treatment Control Unit N2 | Gantry 2 | 400 | 420 | 200 | 15 |
| G303 | Distribution Box TR2 | Gantry 2 | 600 | 800 | 1800 | 450 |
| G307 | Nozzle 2 Electronic Units Cabinet N1 | Gantry 2 | 553 | 780 | 767 | 200 |
| G308 | Nozzle 2 Electronic Units Cabinet N2 (OCE/UA) | Gantry 2 | 553 | 780 | 767 | 200 |
| G309 | Gantry 2 Electronic Cabinet N2 (BPS/BEU) | Gantry 2 | 553 | 780 | 767 | 200 |
| C013 | SICU2 | Pulling Room 2 | 1000 | 320 | 1000 | 70 |
| G306 | Gantry Derivatives Box TR2 - Param Plane1 | Gantry Ph 2 | 1200 | 420 | 1200 | 300 |
| G308 | Gantry Derivatives Box TR2 - Param Plane2 | Gantry Ph 2 | 1200 | 420 | 1200 | 300 |
| G308 | Gantry Derivatives Box TR2 - Control Plane1 | Gantry Ph 2 | 1200 | 420 | 1200 | 300 |
| D032 | Powering Control Unit N3 | Pulling Room 3 | 600 | 800 | 1800 | 450 |
| R351 | X-Ray Generator Cabinet TR3 | Gantry Ph 3 | 500 | 468 | 1101 | 348 |
| G306 | Gantry Driver Cabinet TR3 | Gantry Ph 3 | 800 | 420 | 1000 | F |
| G313 | Robot Controller Unit TR3 (PPS) | Pulling Room 3 | 810 | 600 | 1250 | F |
| C017 | Robot Controller Unit TR3 (PPS) HDC Box | Pulling Room 3 | 600 | 800 | 1800 | 450 |
| D031 | Treatment Control Unit N3 | Gantry 3 | 400 | 420 | 200 | 15 |
| G303 | Distribution Box TR3 | Gantry 3 | 600 | 800 | 1800 | 450 |
| G307 | Nozzle 3 Electronic Units Cabinet N1 | Gantry 3 | 553 | 780 | 767 | 200 |
| G308 | Nozzle 3 Electronic Units Cabinet N2 (OCE/UA) | Gantry 3 | 553 | 780 | 767 | 200 |
| G309 | Gantry 3 Electronic Cabinet N3 (BPS/BEU) | Gantry 3 | 553 | 780 | 767 | 200 |
| C013 | SICU3 | Pulling Room 3 | 1000 | 320 | 1000 | 70 |
| G306 | Gantry Derivatives Box TR3 - Param Plane1 | Gantry Ph 3 | 1200 | 420 | 1200 | 300 |
| G308 | Gantry Derivatives Box TR3 - Param Plane2 | Gantry Ph 3 | 1200 | 420 | 1200 | 300 |
| G308 | Gantry Derivatives Box TR3 - Control Plane1 | Gantry Ph 3 | 1200 | 420 | 1200 | 300 |
| D031 | SICU3 | MCR | 800 | 250 | 600 | 70 |
| L001 | MCR Electronic Cabinet 1 | MCR | 600 | 800 | 1800 | 450 |
| L002 | MCR Electronic Cabinet 2 | MCR | 600 | 800 | 1800 | 450 |
| L003 | MCR Control Cabinet (DAQ) | MCR | 3895 | 6 | 2333 | 1955 |
| A037 | Water Cooling Pump Electric Cabinet | Water Cooling Room | 1500 | 500 | 1400 | 300 |
| A033 | BTS Vault | BTS Vault | 1140 | 800 | 1400 | 660 |
| A034 | Hydraulic Group | BTS Vault | 400 | 800 | 1800 | Floor |



BASEMENT



Remarks :

- Electrical schematics and symbols according IEC standards
- The circuit breakers shall be suitable for isolation and lockable in the "OFF" position. Padlocks to be provided by D/BT.
- D/BT shall provide to PTEV :
 - Complete electrical execution drawings
 - As-built electrical drawings
 - Short-Circuit and Coordination Study of the facility's power distribution system (from city power to hook ups) demonstrating proper coordination of overcurrent protective devices and demonstrating that the short circuit current at each hook up termination is inferior to the corresponding Ampere Interrupting Capacity (AIC) given in the figure above.
 - Arc Flash Study of the facility's power distribution system documenting the available incident energy and the Arc Flash boundary at each hook up circuit breaker as well as at each hook up termination, and documenting the working distance assumed for the calculation.

HK references, see also figure 52.31:

- (A) Hook up 480V Standard
- (B) Hook up 400V Standard (See also table for options on Folio 4/4)
- (C) Hook up 400V Emergency
- (D) Hook up 400V Uninterruptible Emergency

| |
|--|
| PROJECT: PROTON THERAPY |
| S/PROJECT: TATA HBTF MUMBAI |
| V. |
| ELECTRICAL SET |
| TITLE: Main Power distribution Schematic |
| 07.42.33. |
| 52.21 A |

| | | |
|--|-----------------|------------------------------|
| REV: A | DATE: 30/04/15 | MODIFICATION: Original Issue |
| DRAFTSMAN: LCHEN | CHECKED BY: QBA | VALIDATED BY: PV |
| SA. ION BEAM APPLICATIONS | | |
| | | |
| THIS DOCUMENT IS THE PROPERTY OF IBA AND SHALL NOT BE ALTERED COPIED OR COMMUNICATED WITHOUT THE WRITTEN CONSENT OF IBA. | | |

Interface to Control Unit (BLCU or PCU)


The D/BT shall provide an individual interface to the corresponding Control Unit (BLCU or PCU) for each individual main contactor (KM).

Each interface shall include:

- Auxiliary contacts for status information.
- Indicator (e.g., led) which signals the presence of voltage downstream of the main contactor.
- Latch circuit for command of the main contactor. The latch circuit shall be capable to take manual and automatic ON/OFF commands.
- Manual ON/OFF interface (e.g. via momentary push buttons).
- Automatic ON/OFF interface (interface signals shown below).
- Automatic opening of the main contactor (KM) when the upstream main circuit breaker (Q) opens.
- Any command to the main contactor shall be blocked while the main circuit breaker remains open.
- For main contactors connected to X-ray generators, an emergency switch (turn to release type, with red mushroom head) for local emergency OFF action.
- A cable with 12 numbered wires. The cable shall terminate with a loose end (~5m free) to be connected later by the PTEV to the corresponding BLCU or PCU equipment (see figure 52.11). The wires shall be connected as shown below.

Remarks:

- The interface shall be electrically isolated from the rest.
- The main contactor shall not be necessarily located in the same room as the equipment being supplied. But the voltage presence indicator and the manual ON/OFF interface shall be located in the same room as the equipment being supplied.
- For each main contactor connected to an X-ray generator, the voltage presence indicator, the manual ON/OFF interface and the emergency switch shall be located in line of sight of the X-ray generator.
- An example of implementation is given at figure 52.23.



Interface to TSS (SRCU0 and SRCU1)


The D/BT shall provide an individual interface to TSS for each individual main circuit breaker (Q), except for Q8-1.

Each interface shall include:

- Auxiliary contacts for status information.
- Under Voltage Release for emergency OFF.
- A cable with 4 numbered wires. The cable shall terminate with a loose end (~5m free) to be connected later by the PTEV to the corresponding SRCU equipment (see figure 52.11). The wires shall be connected as shown below.

Note:

- The interface shall be electrically isolated from the rest.



| | | | | | | | | | | | |
|------|---|-------|----------|---------------|----------------|------------|-------|-------------|-----|---------------|----|
| REV: | A | DATE: | 30/04/15 | MODIFICATION: | Original Issue | DRAFTSMAN: | LCHEN | CHECKED BY: | QBA | VALIDATED BY: | PV |
|------|---|-------|----------|---------------|----------------|------------|-------|-------------|-----|---------------|----|

MATERIAL: -
SCALE: (A3)
DIMENSIONS: mm
TOLERANCES: -

PROJECT: PROTON THERAPY
SPROJECT: TATA HBTF MUMBAI


V.
ELECTRICAL SET

TITLE:
Interfaces
Requirements

07.42.33.

52.22 A


**EXAMPLE of Interface to
Control System BLCU**



Circuit Breaker equipped with :
 - Auxiliary contacts for status information
 - Under Voltage Release for Emergency OFF

Contactor equipped with :
 - Auxiliary contacts for status information

**EXAMPLE of Interface to
Control System PCU**



Circuit Breaker equipped with :
 - Auxiliary contacts for status information
 - Under Voltage Release for Emergency OFF
 - Ground Fault Protection

Contactor equipped with :
 - Auxiliary contacts for status information

HK for 2nd Xray Gen
 (Reservation for Upgrade/Option as Fluoroscopy or CBCT or Simultaneous xray shots)

| | | | | | |
|--------|----------------|------------------------------|------------------|-----------------|------------------|
| REV: A | DATE: 30/04/15 | MODIFICATION: Original Issue | DRAFTSMAN: LCHEN | CHECKED BY: QBA | VALIDATED BY: PV |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |

MATERIAL: -
 SCALE: (A3) mm
 DIMENSIONS: -
 TOLERANCES: -

PROJECT: PROTON THERAPY
 SPROJECT: TATA HBTF MUMBAI


V.
 ELECTRICAL SET

TITLE:
 Interfaces
 Examples

07.42.33.

52.23 A

THIS DOCUMENT IS THE
 PROPERTY OF IBA AND
 SHALL NOT BE ALTERED
 COPIED OR COMMUNICATED
 WITHOUT THE WRITTEN
 CONSENT OF IBA.



SA. ION BEAM APPLICATIONS

Required Electrical Hookups in the Power Supply Room, BTS and Treatment Rooms

| Hook-up | Load Type | Voltage | Phases | Frequency | Operating Power (a) | Peak Power | Input Current THD (b) | Power Factor | Power Category (c) | INPUT Cables | |
|---------|----------------------------------|------------|---------------------------------|-------------|---------------------|------------|-----------------------|--------------|--------------------|-----------------|--|
| | | Volts | | Hz | kVA | kVA | | | | mm ² | |
| HK1 | Cyclotron Main Coils PS | 400 or 480 | 3 phases + Neutral + Ground (h) | 50 or 60 Hz | 210 | 300(g) | >=15% | >=0,90 | Standard | | |
| HK2 | Cyclotron FPA Anode PS | 400 or 480 | | | 90 | 160(g) | >=10% | >=0,90 | Standard | Standard | |
| HK3 | Dipole PS | 400 | | | 210 | 315(d) | 39% | 0.93 | Standard | Standard | |
| HK4 | Quadrupole PS | 400 | | | 145 | 220 (e) | 6 pulses rectifier | 0.93 | Standard | Standard | |
| HK5 | SMPs | 400 or 480 | | | 80 | N/A | >10% | 0.99 | Standard | Standard | |
| HK6-0 | Cyclo Power Distribution Cabinet | 400 | | | 100 | N/A | | | Standard | Standard | |
| HK7-0 | Cyclo Power Distribution Cabinet | 400 | | | 20 | N/A | | | Emergency | Emergency | |
| HK6-1 | PT Power Distribution Cabinet | 400 | | | 70 | N/A | | | Standard | Standard | |
| HK7-1 | PT Power Distribution Cabinet | 400 | | | 90 | N/A | | | Emergency | Emergency | |
| HK8-1 | PT Power Distribution Cabinet | 400 | | | 60(i) | N/A | | | Uninterruptible | Uninterruptible | |
| HK9 | TR1 x-ray generator | 480 | | | 0,5 (f) | 100 (f) | | | Standard | Standard | |
| HK10 | TR2 x-ray generator | 480 | 0,5 (f) | 100 (f) | | | Standard | Standard | | | |
| HK11 | TR3 x-ray generator | 480 | 0,5 (f) | 100 (f) | | | Standard | Standard | | | |

(a) "Operating Power": This is the power necessary to operate all the equipment at nominal power

(b) THD: Total Harmonic Distortion

(c) "Standard power" is defined as the city power available under normal circumstances. It is used by the system for operating and maintenance purposes.

• "Emergency power" shall always be available, including during a loss of standard power. However, there may be momentary interruption of power service (15 seconds maximum). Emergency power will be used to operate the Gantry and Patient Positioner for purposes of removing the patient following a city power failure.

• "Uninterruptible power" shall always be available, without interruption, even in the event of loss of standard or emergency power. The capacity of the UPS is designed to last 15 minutes at full load. It is required for computers and other critical equipment.

(d) During 0,5s every 60s

(e) During 0,2s every 60s

(f) Each generator (=100KVA) consumes 500 VA stand-by; 120A RMS max. during x-ray shots (less or equal to 1 second). If CBCT or fluoroscopy option up to 15 shots per second.

We propose this table to find gauge wire for X-Ray Hook up cable as a function of length (in accordance with X-Ray Generator Cabinet provider document).

| | |
|----------------------|---|
| Total Length Hook-Up | --> Gauge Cable 3Ph. + N + Ground |
| between 0m to 15m | --> 4x35mm ² + G16mm ² or 4 x AWG2 + Ground |
| between 15m to 30m | --> 4x70mm ² + G35mm ² or 4 x AWG2/0 + Ground |
| between 30m to 45m | --> 4x95mm ² + G50mm ² or 4 x AWG3/0 + Ground |
| between 45m to 60m | --> 4x120mm ² + G70mm ² or 4 x AWG4/0 + Ground |
| between 60m to 75m | --> 4x150mm ² + G95mm ² or 4 x #300MCM + Ground |
| between 75m to 90m | --> 4x185mm ² + G95mm ² or 4 x #350MCM + Ground |

If the D/BT Would like to use another Gauge Cable, the D/BT must produce a calculation sheet which proves that this choice complies with a the total impedance of less than 0.16 Ohm in accordance with IEC 60601-2-7 paragraph 10.2.2.

(g) During more than 60s.

(h) The D/BT shall assure that the 5 wire power cabling shall continue to all IBA loads.

(i) 80 if TPS/OIS system provided by PTEV



S.A. ION BEAM APPLICATIONS

| | | | |
|--------|----------------|------------------------------|------------------|
| REV: A | DATE: 30/04/15 | MODIFICATION: Original Issue | VALIDATED BY: PV |
| | | DRAFTSMAN: -- | CHECKED BY: QBA |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |

MATERIAL: --
 SCALE: (A3)
 DIMENSIONS: mm
 TOLERANCES: --



PROJECT: PROTON THERAPY
 SPROJECT: TATA HBTF MUMBAI

V.
 ELECTRICAL SET

TITLE:
 Required Electrical Hookups Table

07.42.33.

52.24 A



POWER REQUIREMENTS

- The power to be provided is based on the figure 52.24 : Required Electrical Hookups Table D/BT will provide, install and terminate electrical hookups with power and voltage as described in Figure 52.21- Main Power Distribution Schematic at the location shown on this figure
- D/BT shall verify compliance with local regulation and add a quick disconnect box if required.
- All hook-ups terminations shall be at 250 cm($\pm 100''$) from the floor with a loose end ($\sim 5m$ ($\pm 200''$) free) to be connected later by the PTEV to the cabinets.
- Rapid access to the circuit breakers protecting the HK must be provided at all times to PTEV for service purposes. The hookups breakers must disconnect the Neutral conductor line for safety reason.
- The required voltage stability for all the hook-ups is a $\pm 10\%$. Over-voltage suppression will be provided by the D/BT.
- The D/BT shall provide Power OFF button on each x-ray generator hook-up.
 A Power OFF button shall be located in each X-ray control area (near the x-ray console).
 A Power OFF button shall be located in the TCR for PT center having an x-ray console located in TCR (CBCT or fluoroscopy PT center).
 The power OFF button shall power off the X-ray generator hook-up.
 The power OFF button shall be red.
 The power OFF button shall be labeled "X-ray generator hook-up power OFF button"
 The power OFF button shall be connected to the x-ray generator hook-up
 A model of Power OFF button is the XALD112 from Schneider.
- The time delay of the feeder breakers (customer side) shall be adjustable in order to avoid that those breakers trip before the PTEV secondary breakers. The over-current time delay before tripping will be adjustable in the following range: between 15 ms to 230 ms.

THIS DOCUMENT IS THE PROPERTY OF IBA AND SHALL NOT BE ALTERED OR COMMUNICATED WITHOUT THE WRITTEN CONSENT OF IBA.

IBA
ION BEAM APPLICATIONS

| | | | | | |
|-------------|----|---------------|----------|---------------|----------------|
| REV: | A | DATE: | 30/04/15 | MODIFICATION: | Original Issue |
| CHECKED BY: | PV | VALIDATED BY: | | DRAFTSMAN: | LCHEN |
| | | | | CHECKED BY: | OBA |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |

MATERIAL: -
 SCALE: 1/200 (A3)
 DIMENSIONS: mm
 TOLERANCES: -




PROJECT: PROTON THERAPY
 SPROJECT: TATA HBTF MUMBAI

V.
ELECTRICAL SET

TITLE:
 Electrical Hookups Locations

07.42.33.

52.31 A



See also Figure 52.24

UPS REQUIREMENTS

1. UPS outlets for PTEV use only will be located as described on this Figure. These outlets will be installed and wired by D/BT.
2. The UPS outlets will have a color as per local code and a receptacle to prevent accidental connection of non-UPS equipment and to prevent from accidental disconnection of UPS powered equipment.
3. The UPS will be supplied and installed and maintained by D/BT and will comply with IEC62040-1 for safety and IEC62040-2 for EMC. A declaration of conformity from the manufacturer of the UPS certifying the compliance to these 2 standards shall be provided by D/BT.
4. If the UPS is installed in the same room as the PTE, the D/BT should prevent against any risk of contamination of PTE in case of a liquid spillage coming from the batteries.
5. Type: Online - double conversion UPS with static and manual bypass (VFI-SS-111 according IEC 62040-3)
6. Output : 3x400V+N+PE (TN-S)
7. Output rating : see table on figure 52.24
8. Backup time at full load : 15 minutes
D/BT will provide test report proving that the load can be fed for the requested time
9. The cable coming from each individual UPS wall outlet shall terminate with a loose end (≈ 5 m free) to be connected later by PTEV to the equipment power distribution cabinet.
10. Each individual UPS wall outlet shall be on a single circuit
11. CAUTION: extra plugs needed if RT cameras included in configuration

| UPS plugs in other PTEV areas | | | | |
|-------------------------------|----------------|------------|-----------------------------------|--|
| Room | Quantity | Circuit n° | Position | |
| in TPS/OIS server room* | 1x16A | TBD | At the desk On the wall | |
| in TPS/OIS server room* | 4x16A 2x32A | TBD | At the rack Under raised floor | |

* Only if TPS/OIS servers are provided by PTEV

⊗ IEC 60309 2P+E socket outlet with min IP44 protection degree: 16A, except 32A in TPS/OIS server room (see table).
EX: type P17 (Legrand)

THIS DOCUMENT IS THE PROPERTY OF IBA AND SHALL NOT BE ALTERED, COPIED OR COMMUNICATED WITHOUT THE WRITTEN CONSENT OF IBA.

S.A. ION BEAM APPLICATIONS

| | | | | | |
|--------|----------------|------------------------------|------------------|-----------------|------------------|
| REV: A | DATE: 30/04/15 | MODIFICATION: Original Issue | DRAFTSMAN: LCHEN | CHECKED BY: QBA | VALIDATED BY: PV |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |

MATERIAL: -
 SCALE: 1/250 (A3)
 DIMENSIONS: mm
 TOLERANCES: -

PROJECT: PROTON THERAPY
 SPROJECT: TATA HBTF MUMBAI

V.

ELECTRICAL SET

TITLE: UPS connections locations

07.42.33.

52.32 A

| | | | | | |
|-------------|----|---------------|----------|---------------|----------------|
| REV: | A | DATE: | 30/04/15 | MODIFICATION: | Original Issue |
| CHECKED BY: | PV | VALIDATED BY: | | DRAFTSMAN: | LCHEN |
| | | | | CHECKED BY: | OBA |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |

MATERIAL: —
 SCALE: 1/250 (A3)
 DIMENSIONS: mm
 TOLERANCES: —

PROJECT: PROTON THERAPY
 SPROJECT: TATA HBTF MUMBAI

V.
ELECTRICAL SET

TITLE:
 Grounding connections locations




07.42.33.

52.33 A

GROUNDING

1. Grounding (earthing) shall be in accordance with applicable national or local electrical codes.
2. The grounding (earthing) shall be constructed to achieve the lowest practical earth impedance. Multiple ground rods, mesh or grid networks, ground plates, building perimeter ground loops, and other solutions can be used for this purpose. The earth ground resistance shall in any case be less than 10 Ohms.
3. The earth path shall be permanent and continuous. Protective Earth (PE) shall contain no switch or fuse and be such that the disconnection of one assembly does not disconnect other assemblies.
4. Main bonding shall be achieved. Building components, such as structural steel and metallic piping systems, shall be connected to the main electrical ground in accordance with national or local electrical codes.
5. Supplementary bonding shall be achieved. Metallic raceways, cabinet stands, cable trays, metallic water pipes, metallic floor tiles, catwalks, frames, etc. shall be bonded.
6. Grounding busbars shall be installed at the locations described in this Figure to allow for the connection of PTE. Each grounding busbar shall be individually connected to the main earthing terminal with a resistance less than 0.1 Ohm. Each grounding busbar shall also be bonded to the closest accessible building steel.
7. The Neutral (N) and Protective Earth (PE) conductors shall be routed separately as from the transformer or the feeding point and a continuous 5-conductor network L1/L2/L3/N/PE must exist in the building (TN-S electrical distribution network).
8. The D/BT will provide an official inspection report related to the grounding and bonding. The report shall at least contain the measured earth ground resistance, the measured resistance between each grounding busbar and the main earth terminal, and the results of the continuity testing of protective conductors including main and supplementary equipotential bonding.

- W** On the wall 40 cm (1'-4") from finished floor
 - F** On the concrete floor
 - T** In the trench, 20 cm (8") from top of trench, opposite wall of cable tray
 - P** In PPS pit
- Grounding : rectangular copper bar with min. 8 threaded holes M12 (Ø¼")



CABLE TRAYS AND DUCTS

The cables will run from one room to another using ducts or cable trays. This section gives information in order for the D/BT to provide means of getting these cables through the building.

In figure 52.42, a schematic shows the cable distribution in all the Facility and indicates each location where a set of ducts shall begin or end, and references them by a label. It also indicates the quantity and the diameters of the ducts.

All the information provided in this section relates to PTEV needs for PT equipment installation.

Additional needs for equipment not provided by the PTEV, like general lighting, power outlets, CCTV, fire detection ect. are not taken into account.

1. All ducts and cable trays defined hereafter will be for the PTEV use only. No other equipment or cabling can be installed in those without the PTEV written agreement. In that case, ducts or cable trays shall be made larger and their design will be adapted.
2. The detailed design of cable ducts and the routing of cable trays are the responsibility of the D/BT but will be reviewed by the PTEV to ensure compatibility with the cables pulling
3. All ducts and cable trays will be supplied and installed by the D/BT.
4. The Therapy Safety System wiring will be placed in separate and dedicated ducts or cable trays.

CABLE TRAYS REQUIREMENTS

1. All required cable trays must be ladder or wire mesh types, to facilitate getting cables in and out.
2. For all cable trays, wire type is recommended (for example Cablofil CF105 from Cablofil inc.) to minimize installation time and cost. A sufficient number of cable tray supports shall be installed to avoid any bending under the cable weight. If the cable trays are bending due to the cable weight, the D/BT will replace the required parts at its own cost and upon PTEV request. A site survey can be realized by PTEV and the D/BT to list the damaged cable trays.
3. All the cable trays running under a raised floor, will be easily accessible by removing the raised floor panel located just above the tray .
4. The D/BT will avoid sharp edges that could damage the cable insulation.
5. Metallic cable trays shall be grounded. The connection shall be made in such a way that if one section of the tray is removed the other sections are not disconnected. It is recommended to run a ground conductor along the length of the tray, and to connect it to the electrical ground system in at least one point and to each cable tray section via grounding clamps. The applicable national or local electrical codes shall be followed to determine the ground conductor minimum size.
6. All visible cable trays in the MCR, TR and TCR will be closed by the D/BT with an appropriate cover easily removable.

DUCTS REQUIREMENTS

1. Ducts must be designed to minimize radiation leakage.
2. In case the ducts are metal, each duct shall be grounded according to applicable regulations.
3. The D/BT will avoid all cutting edges at the duct ends that could damage the cable insulation.
4. Special care shall also be taken to avoid any obstacles within the tubes. They shall therefore be seamless.
5. All electrical conduits for PTE cabling shall be permanently tagged and numbered at each end by the D/BT, using a convention developed by PTEV. See Figure 52.41 - General comments & legend. The labels will be placed on the wall and on the pulling tapes or cords.
6. All ducts shall be closed by the D/BT at each end during construction to avoid dust or concrete obstructing the ducts.
7. Pulling tapes or cords shall be placed in the ducts from end to end with 100 cm ($\pm 40''$) extra at each end by the D/BT.
8. All ducts passing under the beam line shall be embedded minimum 20 cm ($\pm 8''$) under the concrete surface to avoid interferences with the equipment anchoring. This is also valid for non-PTE ducts.
9. If requested by local regulation, the D/BT will provide fireproof sealing in all PTEV electrical ducts after installation of PTEV cabling and coordinate with PTEV.
10. Electrical conduits ducts end:
 - 10.1 Duct emerging horizontally from the concrete will be minimum at 5 cm ($\pm 2''$) above the ground to avoid water to enter.
 - 10.2 Duct emerging vertically from the floor or ceiling, will exceed at least 5 cm ($\pm 2''$) of the concrete.
11. A list of the conduits lengths or close estimation (± 2 m ($\pm 7'$)) based on figure 52.43 shall be transmitted by D/BT to PTEV 6 months before BOD.

| | | | | | |
|------------|-------|-------------|----------|---------------|----------------|
| REV: | A | DATE: | 30/04/15 | MODIFICATION: | Original Issue |
| DRAFTSMAN: | LOHEN | CHECKED BY: | OBA | VALIDATED BY: | PV |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |

| | |
|-------------|------|
| MATERIAL: | - |
| SCALE: | (A3) |
| DIMENSIONS: | mm |
| TOLERANCES: | - |

PROJECT: PROTON THERAPY
 SPROJECT: TATA HBTF MUMBAI


V. ELECTRICAL SET

TITLE:
 Cable trays & ducts requirements


07.42.33.

52.41 A

LEGEND




Labels of the tubes to be provided by D/BT according to this numbering structure; following a list provided by PTEV:




1. All the Cable trays to be installed by D/BT.
2. The cable trays are dedicated for PTEV equipment only. Other cable trays (type and routing) will be foreseen by D/BT.
3. Coordination and execution drawings will be developed and coordinated with PTEV.


GOOD DUCTS INSTALLATION PRACTICE AND RADIOPROTECTION RECOMMENDATION




Note: Space between 2 ducts must be equal to 2xØ of the ducts to allow the concrete to fill correctly the space in between.




Note: Minimum bending radius of a duct must be equal to 10xØ of duct to ease cable pulling



Note: Maximum 3 bendings along a duct routing



Note: Bind the ducts to steel bars or steel grids into concrete to avoid displacement during the pouring



Note: No or less ducts as possible in an area around the beamline located between 75 cm (±30") to 175 (±70") cm from the floor (50 cm (±20") above and under the beamline). If really needed, after validation by PTEV, the D/BT shall apply strictly the routing and arrangements validated.

Cable Trays

Additional cable trays shall be provided wherever cables will run between 2 locations unlinked by conduits or trenches.

D/BT shall install 1 row of cable tray 100 mm (±4") in the trenches of every treatment room.

To be revised if modification of configuration



| | | | | | |
|------------|-------|-------------|----------|---------------|----------------|
| REV: | A | DATE: | 30/04/15 | MODIFICATION: | Original issue |
| DRAFTSMAN: | LCHEN | CHECKED BY: | QBA | VALIDATED BY: | PV |

| | |
|-------------|------|
| MATERIAL: | --- |
| SCALE: | (A3) |
| DIMENSIONS: | mm |
| TOLERANCES: | --- |

PROJECT: PROTON THERAPY
 SPROJECT: TATA HBTF MUMBAI

V.

ELECTRICAL SET

TITLE:
 General comments & legend

07.42.33.

52.42 A

| | | | | | |
|------------|-------|-------------|----------|---------------|----------------|
| REV: | A | DATE: | 30/04/15 | MODIFICATION: | Original Issue |
| DRAFTSMAN: | LCHEN | CHECKED BY: | QBA | VALIDATED BY: | PV |

MATERIAL: —
 SCALE: (A3)
 DIMENSIONS: mm
 TOLERANCES: —


PROJECT: PROTON THERAPY
 SPROJECT: TATA HBTF MUMBAI

V.
ELECTRICAL SET

TITLE:
 Schematic:
 Cable distribution
 in the building

07.42.33.

52.43 A



----- : Connection made by cable tray

- ⊞ = n x Ø65 mm/Ø2,5" (internal diameter)
- ⊞ = n x Ø100 mm/Ø4" (internal diameter)
- ⊞ = n x Ø25 mm/Ø1" (internal diameter)
- ⊞ = n x Ø200 mm/Ø8" (internal diameter)

To be revised if modification of configuration

See fig. 52.43 for numbering of the sections and of the tubes

THIS DOCUMENT IS THE PROPERTY OF IBA AND SHALL NOT BE ALTERED, COPIED OR COMMUNICATED WITHOUT THE WRITTEN CONSENT OF IBA.

S.A. ION BEAM APPLICATIONS



| | | | | | | | | | | | |
|------|---|-------|----------|---------------|----------------|------------|-------|-------------|-----|---------------|----|
| REV: | A | DATE: | 30/04/15 | MODIFICATION: | Original Issue | DRAFTSMAN: | LCHEN | CHECKED BY: | QBA | VALIDATED BY: | PV |
|------|---|-------|----------|---------------|----------------|------------|-------|-------------|-----|---------------|----|

MATERIAL: —
 SCALE: 1/200 (A3)
 DIMENSIONS: mm
 TOLERANCES: —


PROJECT: PROTON THERAPY
 SPROJECT: TATA HBTF MUMBAI

V.
 ELECTRICAL SET


TITLE:
 Location of cable trays & duct ends
 Treatment Level

07.42.33.

52.44 A



TREATMENT LEVEL




VIEW A

To be revised if modification of configuration


See Fig. 52.43 for numbering of the sections and the tubes

CABLE TRAYS IN PSR & CYCLO PIT

1. The cable trays for the PSR need to be able to support the weight of the cables (40 kg/m (30 lbs/ft)) and workers pulling the cables. Therefore, D/BT will install appropriate cable trays and support framing. See Section A-A for cable trays positioning.
2. In the cyclotron pit, it will consist of 2 cable trays high installed vertically on the wall on 2 levels. See section B-B.




UPPER LEVEL




BASEMENT

SECTION A-A Typical CT Section



SECTION B-B



To be revised if modification of configuration

THIS DOCUMENT IS THE PROPERTY OF IBA AND SHALL NOT BE ALTERED, COPIED OR COMMUNICATED WITHOUT THE WRITTEN CONSENT OF IBA.



S.A. ION BEAM APPLICATIONS

| | | | | | |
|------|----------|----------------|------------|-------------|---------------|
| REV: | DATE: | MODIFICATION: | DRAFTSMAN: | CHECKED BY: | VALIDATED BY: |
| A | 30/04/15 | Original Issue | LCHEN | QBA | PV |

MATERIAL: —
SCALE: 1/200 (A3)
DIMENSIONS: mm
TOLERANCES: —



PROJECT: PROTON THERAPY
SPROJECT:TATA HBTF MUMBAI

V.
ELECTRICAL SET


TITLE:
Location of cable trays & duct ends
Upper Level
Basement

07.42.33.

52.45 A




PULLING ROOM




To be shielded by DBT with high density blocks


VIEW 5




VIEW 6




VIEW 1



VIEW 2



VIEW 3




VIEW 4

| | | | | | |
|--------|----------------|------------------------------|------------------|-----------------|------------------|
| REV: A | DATE: 30/04/15 | MODIFICATION: Original Issue | DRAFTSMAN: LCHEN | CHECKED BY: QBA | VALIDATED BY: PV |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |

MATERIAL: 1/200 (A3)
 SCALE: 1/200 (A3) mm
 DIMENSIONS: mm
 TOLERANCES: -

THIS DOCUMENT IS THE PROPERTY OF LBA AND SHALL NOT BE ALTERED COPIED OR COMMUNICATED WITHOUT THE WRITTEN CONSENT OF LBA.



S.A. ION BEAM APPLICATIONS

PROJECT: PROTON THERAPY
 SPROJECT: TATA HBTF MUMBAI

V. ELECTRICAL SET

TITLE: Cable trays & ducts ends : views and details


07.42.33.

52.46 A


To be revised if modification of configuration

CABLE TRAYS IN MCR


1. The cable trays for the MCR need to be able to support the weight of the cables (20 kg/m, 15 lbs/ft) and workers pulling the cables. For that purpose, D/BT will install appropriate cable trays and support framing. See Section A-A on this figure for cable trays positioning.



Section A-A



THIS DOCUMENT IS THE PROPERTY OF IBA AND SHALL NOT BE ALTERED, COPIED OR COMMUNICATED WITHOUT THE WRITTEN CONSENT OF IBA.



S.A. ION BEAM APPLICATIONS

| | | | | | | | | | | | |
|------|---|-------|----------|---------------|----------------|------------|-------|-------------|-----|---------------|----|
| REV: | A | DATE: | 30/04/15 | MODIFICATION: | Original issue | DRAFTSMAN: | LCHEN | CHECKED BY: | QBA | VALIDATED BY: | PV |
| | | | | | | | | | | | |
| | | | | | | | | | | | |
| | | | | | | | | | | | |
| | | | | | | | | | | | |
| | | | | | | | | | | | |
| | | | | | | | | | | | |
| | | | | | | | | | | | |
| | | | | | | | | | | | |
| | | | | | | | | | | | |

MATERIAL: -
 SCALE: 1/50 (A3)
 DIMENSIONS: mm
 TOLERANCES: -

PROJECT: PROTON THERAPY
 SPROJECT: TATA HBTF MUMBAI

V.
ELECTRICAL SET

TITLE:
 Main Control
 Room cables trays


07.42.33.

52.47 A

CABLE TRAYS IN CYCLO VAULT, ESS/BTS

Considering beam direction as reference, the left-hand side trenches will be dedicated for the electrical service of the beam line.
Plan view and sections of the wiring trenches on the treatment level are shown on this Figure.

1. D/BT shall provide removable walk-on aluminium or galvanized steel metal gratings on all wiring trenches.
Gratings shall be designed not to exceed 25 kg (55 lbs) for manipulation purposes. They will have 10 x 5 cm (4" x 2") cut-outs along grating edges on beam line side for cable and pipe passage from the trench to the equipment (provided by the D/BT). The gratings shall be able to support a moving point load of 1500 kg (3000 lbs).
2. The trenches shall be provided with drains connected to the main floor drain system.
3. Where the trenches must be interrupted due to magnet stands or shielding walls located on their way, the trenches will be connected together by 6 or 8 ducts embedded in the concrete of 100 mm (4") diameter. See cross section A-A on this figure.
4. The cable ducts ending in the trenches will emerge from the side of the trench, exceeding max 2.5 cm (1") the concrete surface.
5. The cross sections view on this figure shows how cable trays shall be installed in the electrical trenches of the Cyclotron/ESS/BTS areas.
6. In the ESS/BTS area, it consists of 3 cable trays of 30 cm (1") width installed on 3 levels (one for control cables, one for signal cables, one for power cables).
7. A free space shall be kept between the cable trays and the beam line side of the trench wall.




SECTION A-A


6 pipes Ø100
in front of the trays

SECTION B-B


3 Trays
300x50




SECTION E-E



SECTION D-D



SECTION F-F



For details of trenches and cyclotron vault pit cover, see figure 01.31-2

| | | | | | |
|------------|-------|-------------|----------|---------------|----------------|
| REV: | A | DATE: | 30/04/15 | MODIFICATION: | Original Issue |
| DRAFTSMAN: | LCHEN | CHECKED BY: | OBA | VALIDATED BY: | PV |

MATERIAL: -
SCALE: 1/150 (A3)
DIMENSIONS: mm
TOLERANCES: -

PROJECT: PROTON THERAPY
SPROJECT: TATA HBTF MUMBAI

V. ELECTRICAL SET

TITLE:
Cable trays in
cyclotron vault,
ESS/BTS

07.42.33.

52.48 A

THIS DOCUMENT IS THE
PROPERTY OF IBA AND
SHALL NOT BE ALTERED
COPIED OR COMMUNICATED
WITHOUT THE WRITTEN
CONSENT OF IBA.
iba
S.A. ION BEAM APPLICATIONS

| | | | | | | | | | | | | | | | | | | | | |
|---------------|----------------|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|
| REV: | | | | | | | | | | | | | | | | | | | | |
| DATE: | 30/04/15 | | | | | | | | | | | | | | | | | | | |
| MODIFICATION: | Original Issue | | | | | | | | | | | | | | | | | | | |
| DRAFTSMAN: | LCHEN | | | | | | | | | | | | | | | | | | | |
| CHECKED BY: | QBA | | | | | | | | | | | | | | | | | | | |
| VALIDATED BY: | PV | | | | | | | | | | | | | | | | | | | |

MATERIAL: ---
 SCALE: 1/150 (A3)
 DIMENSIONS: mm
 TOLERANCES: ---


PROJECT: PROTON THERAPY
 SPROJECT: TATA HBTF MUMBAI

V.
ELECTRICAL SET

TITLE:
 Power Supply
 Room cabinets
 stands


07.42.33.

52.51 A





UPPER LEVEL

View A



View B



| NUMBERING | LENGTH | DEPTH |
|-----------|--------|-------|
| 1 | 5035 | 1100 |
| 2 | 6050 | 850 |
| 3 | 3750 | 450 |
| 4 | 1600 | 850 |
| 5 | 1855 | 950 |
| 6 | 4850 | 650 |

May be divided by the D/B/T

May be divided by the D/B/T

May be divided by the D/B/T

DIMENSIONS IN [mm]

1. Very long stands can be divided in length by D/B/T as indicated next to the table
2. Dimensions for information only, they can be subject to small modifications.
3. All the final dimensions will be confirmed by PTEV minimum 3 months before installation.
4. The closing of the stands above the raised floor is depending on the cooling system selected.

VI. OTHER TECHNICS

60. OTHER TECHNICS

ACRONYMS

| | |
|--------|---------------------------------|
| PT: | Proton Therapy |
| ESS: | Energy Selection System |
| D/BT: | Design Building team |
| PTEV: | Proton Therapy Equipment Vendor |
| BTS: | Beam Transport System |
| TSS: | Therapy Safety System |
| MCR: | Main Control Room |
| SRCU: | Safety Redundant Control Unit |
| TPS: | Therapy planning System |
| OIS: | Oncology information System |
| PSR: | Power Supply Room |
| BOD: | Building Occupancy Date |
| TPS: | Treatment Planning System |
| OIS: | Oncology Information System |
| SRCU: | Safety Redundant Control Unit |
| BTB: | Building Terminal Box |
| NC/NO: | Normally Closed / Normally Open |
| PPS: | Patient Positioning System |

TABLE OF CONTENTS:

| | |
|--|-------------|
| - HVAC Requirements | 61.01 |
| - Sources of Radiations and Shielding | 62.01 |
| - Safety System: Presentation | 63.00 |
| - Safety System: Upper level | 63.01 |
| - Safety System: Basement & Treatment Level | 63.02 |
| - Safety system details | 63.03 |
| - XRay Sign Electrical Wiring | 63.04 |
| - Tables of Operands | 63.05-1/2/3 |
| - Radiation detectors & monitors | 64.01 |
| - Fire detection | 65.01 |
| - Lighting | 66.01 |
| - Network : Presentation | 67.00 |
| - Overall Network topology | 67.01 |
| - Physical Network Topology | 67.02 |
| - Optical fiber connections | 67.03 |
| - Network Requirements Basement | 67.04-1 |
| - Network Requirements Treatment Level & upper Level | 67.04-2 |
| - Communication & Surveillance Equipment | 68.01 |

Refer also to **Chapter II (ROOMS)** for architectural specifications of the rooms



| | | | | | |
|-------------|----|---------------|----------|---------------|----------------|
| REV: | A | DATE: | 30/04/15 | MODIFICATION: | Original issue |
| CHECKED BY: | PV | VALIDATED BY: | | DRAFTSMAN: | LCHEN |
| | | | | CHECKED BY: | QBA |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |

MATERIAL: -
 SCALE: (A3)
 DIMENSIONS: mm
 TOLERANCES: -

PROJECT: PROTON THERAPY
 SPROJECT: TATA HBTF MUMBAI

VI.
OTHER
TECHNICS

TITLE:
 General
 Presentation

07.42.33.

60.00 A